stryker

Endoscopy

45 Liter Insufflator Finding Pneumo

Pneumo Sure[™]





Unparalleled Performance Market Leading Technology

Pneumo Sure™

Increased Accuracy

Stryker's 45L Pneumo Sure provides unparalleled performance, real time pressure sensing and multiple user modes. With its new integrated tube design, the Pneumo Sure manages optimal pressure and flow rates to achieve greater accuracy during procedures. The 45L insufflator is designed to handle a wide variety of cases including pediatrics, bariatrics, and vein harvesting. Its ability to safely maintain pneumoperitenium in a range of conditions brings insufflation in the OR to a new level.



New Touchscreen Interface Quick & Easy Navigation

Interactive user Interface



Software Upgrade

The Pneumo Sure was designed to handle the future needs of the OR through an open architecture software platform.

Touchscreen Interface

The Pneumo Sure's intuitive touchscreen ensures quick and easy navigation through the user menu. The redesigned full color interface was specially designed for surgeons and staff.

Features:

- Higher actual flow rate
- Full color touch panel allows for better viewing angles
- + $\operatorname{SIDNE}^\circ\operatorname{compatibility}$ allows surgeons to control the insufflator
- Change Modes with the touch of a button

Superior Performance Across Multiple Specialties

Exclusive Tubeset Design



High Flow With RTP

- Increased accuracy
- Increased flow
- Ideal for advanced laparoscopic cases



System Specifications

45 Liter PNEUMO SURE™	SYSTEM SPECIFICATIONS
Power Supply Input:	100-240 Vac; 3.15A
Frequency:	50-60 Hz
Power Consumption:	150VA
Dimensions:	5.83" (h) x 18.7" (d) x 12.52"(w)
Weight:	19.84 lbs
Operation Conditions - temp:	50-104 F
Operation Conditions - humidity:	30%-75%

The Pneumo Sure offers 3 different tubesets with the option of using real time pressure for more accuracy and increased flow. With the advancement of RTP, the insufflator is able to continuously pump CO₂ into the abdomen while constantly measuring pressure.

High Flow Tubing

 Increased Flow • Integrated Filter • Ideal for Laparoscopic Cases



Heated With RTP

- Heated CO₂
- Maintains higher flow rates
- Increased accuracy

Product Ordering information

ORDERING INFO.	PRODUCTS
620-040-600	Pneumo Sure™
620-040-610	Pneumo Sure™ XL
620-040-660	High Flow Tubing(disposable)
620-040-680	High Flow Tubing-RTP(disposable)
620-040-690	Heated Tubing-RTP(disposable)
620-030-604	Plume Away 4.0
620-030-606	Plume Away 6.0

stryker

Joint Re	placements	
Trauma,	Extremities & Deform	nities
Craniom	axillofacial	
Spine		
Biologics		
Surgical	Products	
Neuro &	ENT	
Interven	ional Pain	
Navigatio	n	
Endosco	by	
Commun	ications	
Imaging		

EMS Equipment

Stryker Endoscopy 5900 Optical Court San Jose, CA 95138 t: 1 800 624 4422 f: 1 800 729 2917

www.stryker.com

The information presented in this brochure is intended to demonstrate a Stryker product. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with [™] designation are trademarks of Stryker. Products referenced with [®] designation are registered trademarks of Stryker.

Literature Number: **1000900581 (Rev. A)** VID/FWG 2.5k 10/07

Copyright © 2007 Stryker Printed in USA