Feel the advantage.



Cordless Driver 3 Accessories

Sterilization case

 Modular sterilization containers assist in the handling, storage, and sterilization of your Cordless Driver 3 components



Accessories

4100-131-132 Universal Drill Chuck Key

6212-000-000 Small System 6 Battery

6215-000-000 Large System 6 Battery

6110-120-000 System 6 Battery Charger

6126-000-000 System 6 Aseptic Battery Kit

and Battery)

6126-110-000 System 6 Non-sterile Battery

6126-120-000 System 6 Aseptic Housing

6110-120-000 System 6 Battery Charger

6110-625-000 System 6 Battery Module

6102-455-000 System 6 Aseptic Sterilization Case

6126-130-000 System 6 Transfer Shield

(Housing, Transfer Shield,

(6215, 6212, and 6126-110 Batteries)

4200-465-000 Sterilization Case

Ordering Information

Cordless Driver 3

4300-000-000 Cordless Driver 3 Handpiece

Drilling Attachments

4100-062-000 Wire Collet (.028"- .078" diameter) 4100-110-000 Synthes[®] Drill 4100-125-000 Pin Collet (.078" - .125" diameter) 4100-126-000 Adjustable Pin Collet 4100-131-000 1/4" Keyed Drill Chuck 4100-132-000 5/32" Keyed Drill Chuck 4100-133-000 1/4" Keyless Drill Chuck 4100-134-000 1/8" Keyless Drill Chuck 4100-135-000 Hudson[®]/Modified Trinkle Drill 4100-160-000 Trinkle Drill 4100-600-000 1:1 High Speed Bur Attachment 4100-700-000 High Speed Bur Attachment

Reaming Attachments 4100-210-000 Synthes[®] Reamer 4100-231-000 1/4" Keyed Reamer Chuck 4100-235-000 Hudson[®]/Modified Trinkle Reamer 4100-260-000 Trinkle Reamer

Other Attachments

4100-310-000 Right Angle Drive, Synthes® 4100-335-000 Right Angle Drive, Hudson*/Modified 4100-355-000 Radiolucent Right Angle Drive 4100-400-000 Sagittal Saw Attachment 4100-410-000 Radiolucent Adapter

stryker

t Replacements
ma, Extremities, Deformities
iomaxillofacial
e
ogics
ical Products
o & ENT
ventional Pain
gation
oscopy
munications
ing
ent Handling Equipment
Equipment

The information presented in this brochure is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with ™ designation are trademarks of Stryker. Products referenced with ® designation are registered trademarks of Stryker.

Literature Number: 9100-001-124 Rev. None UnDe/P.S.

Copyright © 2007 Stryker Printed in USA

Dat

Stryker Instruments 4100 East Milham Avenue Kalamazoo, MI 49001 USA t: 269 323 7700 f: 800 999 3811 toll free: 800 253 3210

www.stryker.com

stryker

Surgical

CD

The proof is in the performance.



Cordless Driver 3

Battery Charging System

Cordless Driver 3 extends the time-honored tradition of power and versatility in the Sports Medicine arena. With the improved Cordless Driver 3, Stryker continues to offer state-of-the-art battery powered equipment that you have come to expect from Stryker Instruments.



Drilling Attachments




4100-133-000 1/4" Keyless Drill Chuck



4100-134-000 1/8" Keyless Drill Chuck



4100-131-000 1/4" Keyed Drill Chuck




4100-132-000 5/32" Keyed Drill Chuck



4100-700-000 High Speed Bur Attachment



4100-135-000 Hudson® Mod. Trinkle Drill



4100-160-000 Trinkle Drill

4100-125-000


Pin Collet



4100-110-000 Synthes® Drill



4100-600-000 1:1 Bur Attachment



Battery Charger

- Modular design allows for 4 different batteries to charge simultaneously.
- Interactive screen allows user to assess the condition of the battery and read battery usage
- Battery charging station diagnoses battery levels and usage
- Conditioning cycle for discharging non-sterile batteries.

Reaming Attachments



4100-235-000 Hudson® Mod. Trinkle Reamer



4100-260-000 Trinkle Reamer



4100-231-000 1/4" Keyed Reamer Chuck







4100-355-000 Radiolucent Right Angle Drive

4100-210-000 Synthes® Reamer



Battery Options



Aseptic Battery System

- Simple 3 piece aseptic delivery to the sterile field
- 90% more capacity than the System 5 aseptic battery

• Powerful Significant increase in battery capacity

- Simple Easy to handle and deliver in the sterile field
- Versatile Modules accommodate all sterile and non-sterile System 6 batteries

6126-000-000 System 6 Aseptic Battery Kit

Cutting Accessories

 Stryker designs, produces and distributes the industry's highest quality cutting accessories.







4100-310-000 Right Angle Drive Synthes® Reamer



Sagittal Saw Attachment

4100-400-000 Sagittal Attachmer



- Variable speed for precision cutting
- Unique tool-free blade mount reduces blade whip & flex.