Heavy Duty Battery Powered Instruments

stryker

RDLESS DRIVER 2

## **Expanding Your Cordless Capabilities**



A History of Excellence...

# stryker Instruments



# more speed, more torque Driven to Perfection...



More Speed... More Torque...

Proven power and reliability for applications such as:

DC

- Sports Medicine
- Trauma
- Total Joint

#### ...Driven to Perfection

Increased speed and torque are combined with proven cordless reliabitlity to provide unsurpassed power and versatility.



Cordless Driver 2 extends the time-honored tradition of power and versatility in the Sports Medicine arena. With the improved Cordless Driver 2, Stryker continues to offer state-of-the-art battery powered equipment that you have come to expect from Stryker Instruments.



Powerful Capable of tunneling, drilling and reaming.




Versatile


Useful for most orthopaedic procedures. Eliminates the need for multiple handpieces.




**Cordless** Exceptional maneuverability.



50% More Speed 20% More Torque





## Accessories

# Sterilization & Storage/Battery Charging

Modular sterilization container assists in the handling, storage and sterilization of your Cordless Driver 2 components.

Battery charging station diagnoses battery levels and usage.

#### 4200-465 Sterilization Case

• Holds handpiece and all accessories for complete sterilization.

• Convenient inner tray presents instruments on a sterile field.

 Modular design allows for 4 batteries to be charged at the same time.

• Indicator lights provide battery status feedback with "Charge", "Ready", and "Replace" indicators.

> 4110-120 Battery Charger
> Quickly charges Cordless Driver 2 and System 4 batteries.





4100 East Milham Ave Kalamazoo, MI 49001 (800)253-3210 (616)323-7700

For more information, contact your local Stryker Sales Representative.