


Sagittal



5400-100 CORE U Drill



5100-010-334 U Drill Sagittal Saw Attachment



Part Number	Graff Thickness	Max Cut Depth	Description
5100-134-005	5.0mm	20.0mm	Caspar Blade
5100-134-006	6.0mm	20.0mm	Caspar Blade
5100-134-007	7.0mm	20.0mm	Caspar Blade
5100-134-008	8.0mm	20.0mm	Caspar Blade
5100-134-009	9.0mm	20.0mm	Caspar Blade
5100-134-010	10.0mm	20.0mm	Caspar Blade
5100-134-011	11.0mm	20.0mm	Caspar Blade
5100-134-012	12.0mm	20.0mm	Caspar Blade
5100-134-013	13.0mm	20.0mm	Caspar Blade
5100-134-014	14.0mm	20.0mm	Caspar Blade

*Above parts for use with 5100-010-334

Oscillating



5400-31 CORE Oscillating Saw



Part Number	Cut Edge	Thickness	Cut Depth	Description
5400-031-032	11.5mm	.38mm	7.0mm	Precision Thin Intra Oral Blade
5400-031-034	11.5mm	.38mm	12.0mm	Precision Thin Intra Oral Blade
5400-031-050	9.0mm	.38mm	12.0mm	Precision Thin Intra Oral Blade
5400-031-052	5.0mm	.38mm	12.0mm	Precision Thin Intra Oral Blade
5400-031-053	11.5mm	.64mm	7.0mm	Precision Thin Intra Oral Blade
5400-031-054	11.5mm	.64mm	12.0mm	Precision Thin Intra Oral Blade




Part Number	Cut Edge	Thickness	Cut Depth	Description
5400-031-170	9.0mm	.38mm	12.0mm	Crescentic Blade
5400-031-171	13.5mm	.51mm	25.5mm	Crescentic Blade
5400-031-172	18.5mm	.51mm	32.0mm	Crescentic Blade
5400-031-176	13.5mm	.38mm	32.0mm	Crescentic Blade
5400-031-415	9.5mm	.51mm	26.0mm	Crescentic Blade
5400-031-416	15.0mm	.51mm	30.0mm	Crescentic Blade
5400-031-417	18.0mm	.51mm	30.0mm	Crescentic Blade

Reciprocating



5400-37 CORE Reciprocating Saw



Part Number	Thickness	Description
5100-037-901	.38mm	Right IBO
5100-037-902	.38mm	Left IBO