

Reliability for all Laparoscopic Procedures



Laparoscopic Instruments

Reliable Solutions for Pediatric, Bariatric and All Gen




Stryker Endoscopy's laparoscopic instruments are manufactured out of the highest quality stainless steel in conjunction with the most durable modern materials. Each instrument is created using the latest in computer controlled machine technology and perfected by the highest standards of master hand craftsmanship.

3mm



5mm



10mm



General Endoscopic Procedures

3mm

20cm

29cm

5mm

33cm

45cm

10mm

33cm

Milled Radel™ Shaft Insulation Ensures:


- » Durability
- » Patient safety
- » Longer instrument life
- » Reduced maintenance cost

360° Shaft Rotation Provides:

- » One finger control of the instrument
- » Variable resistance based on pressure


Handle Design Facilitates:

- » Better Dexterity
- » Decreased finger fatigue due to larger finger rings



Modular Instrument Design

- » Customization of instruments for each user
- » Easy disassembly for cleaning and maintenance



3mm Laparoscopic Instruments

- Handles
- 20cm Scissors
- 20cm Dissectors

3mm Handles (compatible with both 20cm and 29cm instruments)



PEEK Monopolar Handle
250-282-045



PEEK Multifunction Handle
250-282-046

3mm 20cm Scissors



Curved Metzenbaum Scissors
250-282-000



Peritoneal Scissors
250-282-001



Hook Scissors
250-282-002

3mm 20cm Dissectors



Needle Nose Dissectors
250-282-005



Tapered Bullet Dissectors
250-282-010



Maryland Dissectors
250-282-011

- 20cm Graspers
- 20cm Biopsy
- 20cm Needleholders

3mm

3mm 20cm Graspers



Fundus Grasper
250-282-003



Alligator Grasper
250-282-004



Fenestrated Schertl Grasper
250-282-006



Cobra Grasper
250-282-007



Blunt Grasper
250-282-009



Dorsey Grasper
250-282-096



Double Action Atraumatic Grasper
250-282-100

3mm 20cm Biopsy



Biopsy Forceps
250-282-012



Biopsy Punch
250-282-013

3mm 20cm Needleholders



Straight Needleholder
250-282-014

- Handles
- 29cm Scissors
- 29cm Dissectors

3mm Handles (compatible with both 20cm and 29cm instruments)



PEEK Monopolar Handle
250-282-045



PEEK Multifunction Handle
250-282-046

3mm 29cm Scissors



Curved Metzenbaum Scissors
250-282-047



Peritoneal Scissors
250-282-048



Hook Scissors
250-282-049

3mm 29cm Dissectors



Needle Nose Dissectors
250-282-052



Tapered Bullet
250-282-057



Maryland Dissectors
250-282-058

• 29cm Graspers

3mm

3mm 29cm Graspers



Fundus Grasper
250-282-050



Alligator Grasper
250-282-051



Fenestrated Schertl Grasper
250-282-053



Cobra Grasper
250-282-054



Blunt Grasper
250-282-056



Dorsey Grasper
250-282-097



Single Action Atraumatic Grasper
250-282-099



Double Action Atraumatic Grasper
250-282-101

- 29cm Biopsy
- 29cm Needleholders
- 29cm Monopolar Probes

3mm 29cm Biopsy



Biopsy Forceps
250-282-059



Biopsy Punch
250-282-060

3mm 29cm Needleholders



Straight Needleholder
250-282-061

3mm 29cm Monopolar Electrosurgical Probes



Spatula Tip Probe
250-040-103



L-Tip Probe
250-040-104

5mm Laparoscopic Instruments

- 33cm Handles
- 33cm Scissors

5mm

5mm 33cm Handles



PEEK Monopolar Handle
250-080-617



PEEK Multifunction Handle
250-181-106



Monopolar Handle
250-080-233



Ratcheting Handle
250-080-234

5mm 33cm Scissors



Hook Scissors
250-080-260



Peritoneal Scissors
250-080-262



Straight Metzenbaum Scissors (Long Jaws)
250-080-263



Straight Metzenbaum Scissors
250-080-264



Curved Metzenbaum Scissors (Long Jaws)
250-080-265



Mini Metzenbaum Scissors (Curved)
250-080-266



Endo Metzenbaum Scissors (Curved)
250-080-267



Micro Scissors (Straight)
250-080-269

- 33cm Dissectors
- 33cm Forceps

5mm 33cm Dissectors



Dolphin with Spoon
250-080-280



Tapered Bullet
250-080-281



Maryland Dissector
250-080-282



Curved Kelly Dissector
250-080-283



Tissue Grasper (45°)
250-080-284



Mixer
250-080-285



Needle Nose Dissector
250-080-288



Right Angle Dissector
250-080-694

5mm 33cm Forceps



Meso-Ovarian Forceps
250-080-323



Oviduct Forceps
250-080-324



Tenaculum
250-080-698



Park Forceps
250-080-733

• 33cm Graspers

5mm

5mm 33cm Graspers



Serrated Grasper
250-080-286



Narrow Jaw Grasper
250-080-287



**Single Action
Atraumatic Grasper**
250-080-310



**Double Action
Atraumatic Grasper**
250-080-311



**Fenestrated Schertl
Grasper**
250-080-312



**Blunt Grasper with
Spoon**
250-080-313



DeBakey Grasper
250-080-314



Fundus Grasper
250-080-315



Alligator Grasper
250-080-316



Cobra Grasper
250-080-317



**Serrated,
Fenestrated Grasper**
250-080-318



Dorsey Grasper
250-080-319



**Strong Grasper
(Fenestrated)**
250-080-320



Babcock Grasper
250-080-321



**Grasper, Delicate,
Fenestrated**
250-080-322



Micro Grasper
250-080-325

- 33cm Graspers
- 33cm Biopsy Instrumentation
- 33cm Needleholder Handle and Inserts

5mm 33cm Graspers (continued)



**Allis Grasper
(Single Row of Teeth)**
250-080-591



**Allis Grasper
(Double Row of Teeth)**
250-080-592



Retention Grasper
250-080-697



Endo Grasper
250-080-699



**Paddle Babcock
Grasper**
250-080-700



Fenestrated Grasper
250-080-755

5mm 33cm Biopsy Instrumentation



Biopsy Punch
250-080-290



**Double Action Biopsy
Spoon with Pin**
250-080-291



Blakesly Biopsy
250-080-292



**Single Action Biopsy
Spoon with Pin**
250-080-293

5mm 33cm Needleholder Handle and Inserts



Needleholder Handle
250-080-235



Straight Needleholder
250-080-340



**Curved Left
Needleholder**
250-080-341



**Curved Right
Needleholder**
250-080-342

- Bipolar Handles
- 33cm Bipolar Instrumentation
- 33cm Monopolar Probes

5mm

5mm Bipolar Handles (compatible with both 33cm and 45cm instruments)



Ring Handle
250-080-236




Spring Handle
250-080-237


5mm 33cm Bipolar Instrumentation



Fenestrated Bipolar Forceps
250-080-360



Paddle Bipolar Forceps
250-080-361



Micro Tip Bipolar Forceps
250-080-362



Bipolar Shaft
250-080-371

5mm 33cm Monopolar Electrosurgical Probes



Spatula Tip Probe
250-040-111



L-Tip Probe
250-040-112



J-Tip Probe
250-040-113

- 45cm Handles
- 45cm Scissors
- 45cm Dissectors

5mm 45cm Handles



PEEK Monopolar Handle
250-181-164



PEEK Multifunction Handle
250-181-108



Monopolar Handle
250-080-595



Ratcheting Handle
250-080-596

5mm 45cm Scissors



Hook Scissors
250-080-573



Endo Metzenbaum Scissors (Curved)
250-080-575



Peritoneal Scissors
250-080-751



Curved Metzenbaum (Long Jaws)
250-080-752

5mm 45cm Dissectors



Dolphin with Spoon
250-080-577



Maryland Dissector
250-080-578



Right Angle Dissector
250-080-729



Curved Kelly Dissector
250-080-753

- 45cm Forceps
- 45cm Graspers

5mm

5mm 45cm Forceps



Oviduct Forceps
250-080-587



Park Forceps
250-080-745

5mm 45cm Graspers



Serrated Grasper
250-080-579



**Single Action
Atraumatic Grasper**
250-080-584



**Double Action
Atraumatic Grasper**
250-080-585



Babcock Grasper
250-080-586



Retention Grasper
250-080-725



Endo Grasper
250-080-726



**Paddle Babcock
Grasper**
250-080-727



Fenestrated Grasper
250-080-732

- 45cm Graspers
- 45cm Biopsy Instrumentation
- 45cm Needleholder Handle and Inserts

5mm 45cm Graspers (continued)



Fundus Grasper
250-080-735



Dorsey Grasper
250-080-736



DeBakey Grasper
250-080-748



**Allis Grasper
(double row of teeth)**
250-080-754

5mm 45cm Biopsy Instrumentation



Biopsy Punch
250-080-580



**Double Action Biopsy
Spoon with Pin**
250-080-581

5mm 45cm Needleholder Handle and Inserts



Needleholder Handle
250-080-235



Straight Needleholder
250-080-588



**Curved Left
Needleholder**
250-080-589



**Curved Right
Needleholder**
250-080-590

- Bipolar Handles
- 45cm Biopsy Instrumentation
- 45cm Monopolar Probes

5mm

5mm Bipolar Handles (compatible with both 33cm and 45cm instruments)



Ring Handle
250-080-236




Spring Handle
250-080-237

5mm 45cm Bipolar Instrumentation



Fenestrated Bipolar Forceps
250-080-350



Paddle Bipolar Forceps
250-080-351



Bipolar Shaft
250-080-370

5mm 45cm Monopolar Electrosurgical Probes



Spatula Tip Probe
250-040-116



L-Tip Probe
250-040-117



J-Tip Probe
250-040-118

10mm Laparoscopic Instruments

- 10mm, 33cm Handles
- 10mm, 33cm Scissors
- 10mm, 33cm Graspers, Forceps, and Dissectors

10mm, 33cm Handles



PEEK Monopolar Handle
250-080-618



PEEK Multifunction Handle
250-181-107



Monopolar Handle
250-080-231



Ratcheting Handle
250-080-230

10mm, 33cm Scissors



Operating Scissors, Serrated
250-080-247

- 10mm, 33cm Graspers, Forceps, and Dissectors

10mm

10mm, 33cm Graspers, Forceps and Dissectors



Allis Grasper
250-080-240



Babcock Grasper
250-080-241



Pennington Grasper
250-080-242



Koeln Grasper
250-080-243



Single Action Claw
250-080-244



Double Action Claw
250-080-245



Biopsy Spoon
250-080-246



Maryland Dissector
250-080-248




Mixer
250-080-249



Paddle Babcock Grasper
250-080-250

Sterilization Trays

- Laparoscopic Instrument Tray
- Metal Laparoscopic Instrument Tray
- 45cm Instrument Tray



Laparoscopic Instrument Tray

Compatible with 3mm, 5mm, and 10mm instruments
250-015-300



Metal Laparoscopic Instrument Tray

250-015-600



45cm Instrument Tray

Compatible with 5mm and 10mm instruments
250-015-704



