

# Versatility

At its Finest

**Mainstay Zip Anchor**



# Mainstay Zip Anchor

Versatility At Its Finest.

## Convenience

---

The Mainstay Zip Anchor consists of a titanium anchor, suture, and needles pre-assembled on a sterile single-use inserter. Three anchor sizes are available: 2.7mm, 3.5mm, and 4.5mm.



4.5mm Mainstay Zip Anchor



4.5mm Stop Drill

## Versatility

---

Mainstay Zip Anchors are pre-assembled. The pre-assembled anchor is convenient for the surgeon and staff and requires no set-up. The Zip anchor can also be repositioned.



3.5mm Mainstay Zip Anchor



3.5mm Stop Drill

## Choices

---

The anchor, suture, and cutting needles are available on a single-use inserter. Zip Anchors come pre-assembled with two strands of braided polyester suture. To facilitate ease of knot tying, one strand of suture is colored green and one strand is colored white.



2.7mm Mainstay Zip Anchor



2.7mm Stop Drill

## An Exceptional Anchor

With each Mainstay Zip Anchor, you get an excellent implant design. The “cradle” is designed to maximize surface area contact and stress distribution potentially increasing suture strength. High pull-out strength, increased suture integrity, and removability make the Zip Anchor an exceptional implant.



## Ordering Information

2.7mm Zip Anchor with #0 suture  
**3910-002-027**

3.5mm Zip Anchor with #2 suture  
**3910-002-035**

3.5mm Zip Anchor with two strands of #0 suture  
**3910-002-235**

4.5mm Zip Anchor with two strands of #2 suture  
**3910-002-245**

4.5mm Zip Anchor with #5 suture  
**3910-002-545**

Soft Tissue Protector - Slotted  
**3910-003-010**

1.5mm Stop Drill (2.7mm anchor)  
**3910-004-015**

2.0mm Stop Drill (3.5mm anchor)  
**3910-004-020**

2.4mm Stop Drill (4.5mm anchor)  
**3910-004-024**

**All suture is pre-loaded with cutting needles.**

## Indications

The Mainstay Soft Tissue Anchor is intended only for the fixation of non-absorbable braided polyester surgical suture to bone for the anatomical locations listed here:

### Shoulder

- Bankart repair
- SLAP lesion repair
- Biceps tenodesis
- Acromio-clavicular separation repair
- Rotator cuff repair
- Capsular shift/capsulolabral reconstruction
- Deltoid repair

### Pelvis

- Bladder neck suspension procedures

### Foot and Ankle

- Medial or lateral instability repair/reconstruction
- Achilles tendon repair/reconstruction
- Mid-foot reconstruction
- Hallux valgus reconstruction

### Elbow, Wrist, and Hand

- Scapholunate ligament reconstruction
- Ulnar/radial collateral ligament reconstruction
- Tennis elbow repair
- Biceps tendon re-attachment

### Knee

- Extra capsular repairs:
  - Medial collateral ligament reconstruction;
  - Lateral collateral ligament;
  - Posterior oblique ligament;
  - Illiotalband tenodesis
- Patellar tendon repair

### Hip

- Abductor tendon reconstruction



---

**Joint Replacements**

---

**Trauma, Extremities & Deformities**

---

**Craniomaxillofacial**

---

**Spine**

---

**Biologics**

---

**Surgical Products**

---

**Neuro & ENT**

---

**Interventional Pain**

---

**Navigation**

---

**Endoscopy**

---

**Communications**

---

**Imaging**

---

**Patient Handling Equipment**

---

**EMS Equipment**

325 Corporate Drive  
Mahwah, NJ 07430  
t: 201 831 5000

[www.stryker.com](http://www.stryker.com)

A surgeon should always rely on his or her own professional clinical judgment when deciding to use which products and/or techniques on individual patients. Stryker is not dispensing medical advice and recommends that surgeons be trained in orthopaedic surgeries before performing any surgeries.

The information presented is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Mainstay and Stryker. All other trademarks are trademarks of their respective owners or holders.

Literature Number: 1000-900-316C  
MS/GS 03/08

Copyright © 2008 Stryker  
Printed in USA