

SwitchPoint Infinity 3

Orchestrating the OR



Simplified

OR Control

Stryker's SwitchPoint Infinity 3: the central point for managing your OR video and surgical devices.

The intuitive interface allows you to easily control your Operating Room through features like one touch routing and video multi-viewing capabilities.



The OR Integration System

Tailored to Your Needs

- Any-to-Any Digital Video Signal Routing
- Highly Configurable and Extensible
- Designed Specifically for Use in the OR
- Built For Future Expansion



Supported Video Standards

- S-Video
- Composite
- RGBHV
- 3G-SDI
- HD-SDI
- DVI

The SwitchPoint Infinity 3 is compatible with standard video signals and supports resolutions including 1080P and up to 1200P.

Your Control Center in the Operating Room

Designed specifically for the operating room, the SwitchPoint Infinity 3 features tools that help staff perform tasks efficiently.

Surgical Checklist Integration

Engage team members during surgical time-outs by broadcasting a pre-surgical checklist to all monitor screens.

Customizable Room Presets

Save and preview frequently used room configurations to reduce turn-around time and expedite room setup.

Split Screen and Quad View

Display up to four video sources for viewing multiple modalities on a single monitor.

Centralized Suite Control

Control compatible Stryker surgical equipment from the touch panel to perform tasks and assist team members from a single location.

Control Multiple Computers

Access and control OR computers using a single keyboard and mouse, saving space and reducing clutter.

Teleconference from the OR

Bridge the distance between sites across campus or locations abroad using IP and baseband video-conferencing features.



World-Class Service and Support

Powered by SORN, Stryker's comprehensive system for remote device management, the SwitchPoint Infinity 3 automatically sends system alerts and maintenance information to technical staff and the Stryker Customer Care team.

SORN immediately notifies Stryker if issues arise, allowing us to respond proactively and keep the system ready to operate when you are.

Remotely Monitor

- Software status
- Hardware maintenance
- Temperature and fan status
- Software versions

Joint Replacements

Trauma, Extremities & Deformities

Craniomaxillofacial

Spine

Biologics

Surgical Products

Neuro & ENT

Interventional Spine

Navigation

Endoscopy

Communications

Patient Care & Handling Equipment

EMS Equipment

Stryker Communications
1410 Lakeside Parkway #100
Flower Mound, TX 75028 (USA)
Tel: (866) 841-5663

www.stryker.com

The information presented in this brochure is intended to demonstrate a Stryker product. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: SORN, Stryker, SwitchPoint Infinity. All other trademarks are trademarks of their respective owners or holders.

Literature Number: 1000900669 (Rev. C)
BC/EX 2000 12/11

Copyright © 2011 Stryker
Printed in USA

www.stryker.com/spi3

