# CORE Sumex Drill

Your Specialty Is Our Priority.

CO Somet Date

Combining the best of our power and ergonomic designs to give you our best, high performance drill.

### **stryker**





MIS Burs

#### Zyphr Burs

- Free cutting; efficient, effective bone removal with less effort
- Smoother cuts and reduced chatter
- Flute technology maximizes bone chip clearance and reduces clogging
- · Neck design improves line of sight



#### **Sumex Drill**

#### 5400-130-000

- 75000 RPM
- Electronic Torque Feedback (Increases RPM's when the drill is engaged in more demanding tasks)
- Tapered front end (Allows for better line of sight)
- Takes SD/PD Attachments and Maestro Attachments
- Fixed and Steering Duragaurds
- · Zyphr, Elite and MIS burs
- Complete OR Standardization





#### **Ordering Information**

<b>Attachments</b>	
5100-120-450	Saber Straight Attachment
5100-120-452	Saber Angled Attachment
5100-120-470	Saber Long Straight Attachment
5100-120-472	Saber Long Angled Attachment
5100-120-480	Saber Extra Long Straight Attachment
5100-120-482	Saber Extra Long Angled Attachment
5100-120-490	Saber MIS Attachment
5100-120-492	Saber Angled MIS Attachment
5400-210-020	Maestro Straight O Attachment
5400-210-050	Maestro Straight M Attachment
5400-210-070	Maestro Straight AM Attachment
5400-210-080	Maestro R Attachment

5400-210-090	Maestro Rx Attachment	
5400-210-060	Maestro Perforator Chuck	
5400-210-501	Maestro Adjustable Depth Drill Guide	
5400-210-011	Maestro Right Angle Attachment	
MIS Attachments		
5100-120-920	SD/PD MIS Medium Straight	

0100 120 720	OD/I D IIIIO IIICaiaiii otiaigiit
5100-120-950	SD/PD MIS Long Straight
5100-120-922	SD/PD MIS Medium Curved
5100-120-952	SD/PD MIS Long Curved
Duraguards:	
5400 210 057	Magetro Small Fixed Duraguar

Duraguards:	
5400-210-057	Maestro Small Fixed Duraguard
5400-210-058	Maestro Medium Fixed Duraguard
5400-210-059	Maestro Large Fixed Duraguard

5400-210-257	Maestro Small Steering Duraguard
5400-210-258	Maestro Medium Steering Duraguard
5400-210-259	Maestro Large Steering Duraguard
5400-210-056	Maestro D-Attachment

## 5400-071-057 Small Duraguard Router 5400-071-058 Medium Duraguard Router 5400-071-059 Large Duraguard Router

Accessories:	
5400-050-000	Core Console
5100-008-000	TPS Footswitch
5400-131-000	Hand Switch

**Routers:** 

#### This document is intended solely for the use of healthcare professionals.

The information presented in this brochure is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with  $^{\text{\tiny TM}}$  designation are trademarks of Stryker. Products referenced with  $^{\text{\tiny B}}$  designation are registered trademarks of Stryker.



Literature Number: MTX9100001060 MTX/GS 06/07 Stryker SA Cité-Centre Grand-Rue 90 1820 Montreux Switzerland

t: +41 21 966 12 01 f: +41 21 966 12 00

www.europe.stryker.com