

# System 8 power tools



# A family of power tools **built for the whole hospital**

# Stryker's System 8 The tools you trust – now more reliable than ever

Stryker's System 8 power tools deliver more of the long-lasting reliability you trust. Built for your O.R. Built for your sterile processing department. And built to benefit your bottom line. In other words, built for the whole hospital.

So put each power tool through its paces. We did, so that we know when the System 8 set passes from our hands to yours, you can have **Total Confidence**.



## Total Confidence



Sports medicine and ortho trauma

Orthopaedic and cardiothoracic surgery



# Meet Stryker's

# System 8 power tools











cordless driver gives you the confidence you need to drill, cut and ream in sports med and ortho trauma procedures. Plus, our exclusive PoweReam Technology gives you a significant increase in torque, for tackling challenging reaming procedures.

System 8 SABO Sagittal Saw – The next iteration of Stryker's SABO Saw, the lightweight System 8 SABO Sagittal Saw delivers optimized cutting speed, control and accuracy in every case.

Stryker Precision Saw – The oscillating tip saw that reduces noise and vibration, enhancing handpiece control for a more precise surgical experience.

#### An award-winning design, made even better

Deliberate, substantive changes give System 8 handpieces improved ergonomics, and give you a better wrist position and firmer grip.

#### Consistent power and uptime

High-capacity cells deliver improved performance to get you through the toughest cases. Plus, large batteries now last more than 2 times longer over the course of their total lifespan without needing replacement.

#### A keyless chuck as convenient as it is reliable

Quick and efficient, the System 8 locking keyless chuck has been engineered to prevent loosening, with the addition of a secondary locking mechanism.

Stryker's EZout Powered Acetabular Revision System is the power tool system that gives you precise, clean cuts

### A pin collet built to last

Advanced material and coating helps to prevent stick and slip. Plus, they're now designed for both threaded and fluted pins, as well as straight.

### A sag saw that boosts uptime in the O.R.

Ouality improvements to the System 8 sagittal saw greatly increase reliability, minimizing product downtime.

### An oscillating tip saw that will evolve the way you cut

Our latest evolution of the Stryker Precision Saw features an enhanced design and expanded line of Precision Falcon Cartridges for ankle, shoulder, hip and knee arthroplasty.

### Long-lasting batteries with increased heat resistance

System 8 batteries can live on the shelf for 30 days without compromising performance. Plus, insulation within each battery means increased resistance to thermal damage during extended exposure to high temperatures.

#### More water resistant than ever

System 8 handpieces are built to be actively washed – able to be temporarily submerged when you need to do so during cleaning, prior to sterilization.

#### **Optimize your power tool management**

Originated through research and guided by our customers, we are pleased to offer Smart Equipment Management (SEM), a revolutionary, online data module designed to help you optimize the management of your System 8 power tools through device-derived, actionable data.

- Proactive device maintenance
- Balanced product use
- Spend justification
- Simplified service
- Asset tracking















## System 8 power tools

# Built for the **whole**

For questions or to place your order today, please contact your Surgical sales representative or call 800 253 3210.

#### Stryker's System 8 part numbers

#### Handpieces

	8202-999-000	EZout System	30	
	8202-000-000	EZout Saw		
	8203-000-000	Single trigger rotary drill	30	
	8205-000-000	Dual trigger rotary drill	30	
	8206-000-000	Reciprocating saw		
	8207-000-000	Sternum saw	30	
	8208-000-000	Sagittal saw		
	8209-000-000	Stryker Precision Saw	30	
	4505-000-000	Cordless driver	30	
	4508-000-000	SABO Sagittal Saw	0	
	Batteries			
	7110-120-000	In increal bettern aborren	C	
	8212-000-000	Universal battery charger System 8 battery pack, small	4	
	8212-000-000	System 8 battery pack, large		
	0213-000-000	System o battery pack, large	4	
Smart Equipment Management 4				
	0280-100-010	SEM user license fee (qty 1)	4	
	0280-200-100	SEM handpieces - data subscription	44 44	
	0280-200-105	SEM batteries - data subscription	44	
	0280-200-120	SEM universal charger - data subscription	44	
		-	44	
Automatic cleaners and detergents 4				
	3000-009-000	Auto62 Enzymatic Cleaner; non-foaming high-rinsing	41	
		enzymatic concentrate	41	
		2/2.5 gallon (9.46L)	41	
	3000-010-000	Auto62 Neutral Detergent; low-foaming high-rinsing	41	
		neutral concentrate	41	
		2/2.5 gallon (9.46L)	4	
	3000-011-000	Auto62 Instrument Lubricant; neutral non-silicone	4	
		instrument lubricant	41	
		2/2.5 gallon (9.46L)	4	

#### **Manual cleaners and detergents**

3000-003-000	Blu62 Pretreatment Foam - 10/32 oz (947mL)
3000-007-000	Blu62 Pretreatment Foam refill - 4/1 gallon (3.78L)
3000-002-000	Blu62 Instrument Detergent - 20/8 oz (237mL)
3000-004-000	Blu62 Instrument Detergent - 3/100 oz (2.97L)
3000-005-000	Blu62 Instrument Detergent - 2/2.5 gallon (9.46L)
3000-001-000	Blu62 Instrument Detergent Dispenser

#### Cordless driver attachments

4505-133-000	1/4" (6.4mm) Locking keyless chuck
4505-134-000	5/32" (4.0mm) Locking
4505-154-000	keyless chuck
4100-131-000	1/4" Drill chuck with key
1100 101 000	•
4100-132-000	5/32" Drill chuck with key
4405-210-000	PoweReam large AO
4405-213-000	PoweReam Hudson
4405-231-000	PoweReam 1/4" chuck with key
4405-235-000	PoweReam Hudson
	modified trinkle
4405-260-000	PoweReam trinkle
4100-126-000	Adjustable pin collet
	2.0 - 3.2mm
4100-062-000	Wire collet 0.7 - 1.8mm
4100-125-000	Pin collet 2.0 - 3.2mm
4100-225-000	Pin collet 3.0 - 4.2mm
4100-110-000	Small AO drill
4100-120-000	T-Latch
4100-135-000	Hudson modified trinkle drill
4100-160-000	Trinkle drill
4100-210-000	Large AO reamer
4100-231-000	l/4" Reamer chuck with key
4100-235-000	Hudson modified trinkle reamer

4100-260-000 Trinkle reamer 4100-215-000 DHS/DCS reamer 4100-355-000 Radiolucent right angle drive 4100-400-000 Sagittal saw attachment 4100-410-000 Right angle drive adaptor 4100-600-000 1:1 Bur attachment 4100-700-000 High speed bur attachment 4100-131-132 Universal chuck key 1331-001-009 1/4" Chuck key 1355-005-000 5/32" Chuck key

#### **Rotary attachments**

8203-026-000	Dual trigger wire collet 0.7 - 2.0mm
8203-036-000	Single trigger wire collet 0.7 - 2.0mm
8203-126-000	Dual trigger pin collet 2.0 - 3.2mm
8203-136-000	Single trigger pin collet 2.0 - 3.2mm
8203-226-000	Dual trigger large collet 3.0 – 4.2mm
8203-236-000	Single trigger large collet 3.0 – 4.2mm
8203-133-000	l/4" (6.4mm) Locking keyless chuck
6203-110-000	Small AO quick connect
6203-113-000	Hudson quick connect
6203-131-000	l/4" Chuck with key
6203-132-000	5/32" Chuck with key
6203-135-000	Hudson modified trinkle quick connect
6203-150-000	Modified trinkle quick connect
6203-160-000	Trinkle quick connect
6203-210-000	Large AO quick connect
6203-215-000	DHS/DCS quick connect
6203-170-000	Triathalon adaptor
7207-003-000	Sternum blade guard
7207-002-000	Sternum blade guard long

#### **ProCare Services**

Our ProCare Services solutions help optimize the reliability and life of your equipment. Our service plans are designed to make it easier for you to get the most out of your investment in equipment from Stryker with options ranging from unlimited priority repairs to preventive maintenance. Please contact your Surgical sales representative for detailed ProCare Services options.

#### **Cleaners and detergents**

Protect your investments with our cleaning and detergent solutions. **Blu62 Instrument Detergent** – a cleaner that protects instruments from the corrosive effects of water and water impurities, as well as corrosion from mixed metal. **Blu62 Pretreatment Foam** – a spray formulated to start the cleaning process in the O.R. and prevent the drying of residue on soiled instruments. **Auto62 Cleaning Chemistries** – an enzymatic cleaner, neutral detergent and instrument lubricant – exclusively designed for use in the automatic wash cycle to remove tough soils with consistent and effective results.

The information presented is intended to demonstrate the breadth of Stryker's product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker's products in your area.

🔟 See Stryker's System 8 power tools Instructions for Use (IFU) for complete details.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Auto62, Blu62, EZout, PoweReam, Precision Falcon, ProCare, SABO, Stryker, Stryker Precision, Total Confidence. All other trademarks are trademarks of their respective owners or holders.

© 2017 Stryker | www.stryker.com

Stryker Instruments 4100 East Milham Avenue Kalamazoo, MI 49001 USA t: 269 323 7700 f: 800 999 3811 toll free: 800 253 3210

Literature number: 9100-003-973 Rev none PF/PS