Heavy Duty Battery Powered Instruments

Modular sterilization container assists in the handling, storage, and sterilization of your System 5 components.

#### 4102-450 Sterilization Case

- Holds 3 handpieces and all accessories for complete
- Convenient inner tray presents instruments on a sterile field.

### 4110-120 Battery Charger

- Quickly charges System 5 and Cordless Driver 2 batteries.
- Non-sterile battery (red) available
- Conditioning cycle completely discharges non-sterile battery before charging.
- Modular design allows for 4 different batteries to be charged at the same time.
- Indicator lights provide battery status feedback with "Charge", "Ready", and "Replace" indicators.



Battery charging station diagnoses battery levels

and usage.



# Complete System 5 Ordering Information



### System 5 Handpieces

PART #	DESCRIPTION
4203	Rotary Handpiece, Single Trigger
4205	Rotary Handpiece, Juligier Higger Rotary Handpiece, Dual Trigger
4206	Reciprocating Saw
4208	Sagittal Saw
Drilling	Attachments

# Drilling Attachments

4103-82-131	Keyless Jacobs® Chuck Adapter
4103-110	Synthes® Chuck
4103-126	Pin Collet
4103-131	1/4" Drill with Jacobs® Chuck (1:1)
4103-133	Keyless Jacobs® Drill Attachment (1:1)
4103-135	Hudson®/Modified Trinkle Drill
4103-160	Trinkle Drill Attachment (1:1)
4103-180	High Torque 1/4" Drill with Jacobs® Chuck (2:1)

# Reaming Attachments

PART #	DESCRIPTION
4103-210 4103-213	Synthes® Reamer (3.25:1) Hudson® Reamer (3.25:1)
4103-231	1/4" Reamer with Jacobs® Chuck (3.25:1)
4103-235 4103-260	Hudson®/Modified Trinkle Reamer (3.25:1) Standard Trinkle Reamer (3.25:1)
4103-413 4103-435	High Torque Hudson® Adapter High Torque Hudson®/Modified Trinkle Reamer (4:1)


# System 5 Accessories

PART #	DESCRIPTION
4115	System 5 Battery
4126	System 5 Aseptic Battery Kit
	(includes battery pack, housing, & shield)
4126-110	Non-Sterile Large Battery
4126-120	System 5 Aseptic Battery Housing
4126-130	Large Transfer Shield 4110-115 Battery Adapter
	(for use with 4115 batteries in System 2000 charger)
4110-120	Battery Charger, 120 Volt
4110-215	Charger Module (2115 battery in System 5 Charger)
4110-415	Charger Module (4115 battery)
4110-426	Non-Sterile Battery Charger Module (4126-110 battery)
5100-901-906	Electrical Connection Brush
2106-151	Sterilization Case (holds 4 handpieces)
2106-152	Sterilization Basket (holds 2 handpieces)
4102-450	Plastic Sterilization Case (holds 4 handpieces)
4102-451	Sterilization Case (holds 3 handpieces)
4102-452	Mini Sterilization Case (holds 2 handpieces)

## **Sternum Applications**

	• • • • • • • • • • • • • • • • • • • •
PART #	DESCRIPTION
4207	Reciprocating/Sternum Saw
4112	Cordless Driver Battery (for use with 4207 Sternum Saw)
4212	Small Smart Cordless Driver Battery (for use with 4207)
4215	Large Smart Cordless Driver Battery (for use with 4207)
4110-412	Charger Module (4112 battery)
4110-422	Cordless Driver Smart Charger Module
	(4212 & 4215 batteries)
4107-008	Sternum Blade Guard
298-97-100	Sternum Blade*
298-97-101	Aggressive Sternum Blade





- Non-sterile battery ensures longer life of battery.
- Battery is easily charged in System 5 charger.
- Simple assembly for operating room







\* See Bur and Blade literature for complete listing of available cutting accessories

# **stryker**° **INSTRUMENTS**

For more information, contact your local Stryker Sales Representative.

4100 East Milham Ave Kalamazoo, MI 49001 (800)253-3210 (269)323-7700

stryker° INSTRUMENTS System 5 Power of Choice
dual trigger, dual speed
for greater versatility

# System 5 Rotary Handpieces

Stryker Instruments, the leader in cordless power, brings its 5th generation of battery powered instrumentation...



#### Collets **Drilling & Reaming Attachments**

1:1 1/4" Jacob's® Drill

Hudson® Modified Trinkle Drill

4103-160 1:1 Standard Trinkle Drill

4103-133 Keyless Drill

4103-82-131 Keyless Chuck Adapter

4103-110 1:1 Synthes® Drill

4103-131

4103-135

Colored coded rings on drilling and reaming attachments allow for easy recognition and selection.


Stryker Instruments offers a wide assortment of quick connect attachments to be used with a variety of reamers, drill bits, and cutting accessories.



4203 Rotary Handpiece

Rotary handpieces allow for both high speed drilling and high torque reaming.

quick and easy selection of direction/safety options.



# System 5 Heavy Duty Saws

Advanced Motor Technology provides increased torque for a faster and more efficient cut.



