

stryker®

Endoscopy

Durable Precision

At Your Fingertips

Conquest™ Arthroscopic Manual
Instruments



Conquest™

Durable Precision

Stryker Endoscopy brings to market the premium line of arthroscopic manual instruments: Conquest™. Designed for a broad range of applications, Conquest™ instruments meld world-class craftsmanship with unparalleled technical precision. Stryker offers a complete line of 40 instruments including Punches, Scissors, Graspers, and Small Joint Instruments.

Manufactured to the tightest tolerances of Electrical Discharge Machining, Conquest™ tips cut with precision, case after case.



Milled from a solid piece of hardened steel and coated in chromium for increased protection, Conquest™ moving parts boast exceptional strength yet provide fluid motion when cutting or biting.

Ergonomically designed to minimize fatigue, Conquest™ handles cut with greater ease and power.

Conquest™ Punches

These instruments are available in a variety of tip configurations to accommodate varying needs. Stryker engineers have designed these instruments with teeth that allow the surgeon to smoothly pull tissue into the jaws of the instrument, providing safe and precise resection.



3.4mm Backbiter
300-034-800



3.4mm Straight Big Bite Punch
300-034-100



3.4mm 30° Right Big Bite Punch
300-034-101



3.4mm 30° Left Big Bite Punch
300-034-102



3.4mm 15° Up Big Bite Punch
300-034-103

2.7mm 30° Right Big Bite Punch
300-027-101

2.7mm 30° Left Big Bite Punch
300-027-102

2.7mm 15° Up Big Bite Punch
300-027-103



3.4mm Upshaft Big Bite Punch
300-034-104



3.4mm 45° Left Tip Biter
300-034-110



3.4mm Straight Scoop Jaw Punch
300-034-300

2.7mm Upshaft Big Bite Punch
300-027-104


2.7mm 45° Right Tip Biter
300-027-111



3.4mm 30° Right Scoop Jaw Punch
300-034-301



3.4mm Narrow Punch
300-034-200



3.4mm 15° Up Narrow Punch
300-034-203

3.4mm 30° Left Scoop Jaw Punch
300-034-302

2.7mm Narrow Punch
300-027-200



3.4mm Right Meniscal Crescent Punch
300-034-500



3.4mm 30° Left Meniscal Crescent Punch
300-034-502



3.4mm 90° Right Rotary Punch
300-034-400

3.4mm Left Meniscal Crescent Punch
300-034-501

3.4mm 30° Right Meniscal Crescent Punch
300-034-503

3.4mm 90° Left Rotary Punch
300-034-401

Conquest™ Scissors

The smooth nose of these scissors allows for safe entry and reduces the likelihood of damaging tissue.



3.4mm Hook Scissors
300-034-600



3.4mm 30° Left Hook Scissors
300-034-602



3.4mm 30° Right Hook Scissors
300-034-601



3.4mm Micro Scissors Punch
300-034-603

2.7mm Micro Scissors Punch
300-027-603

Conquest™ Small Joint Instruments

The 2.0mm diameter makes these instruments the ultimate tool for elbow, wrist, and ankle arthroscopy.



2.0mm Small Joint Punch
242-010-200




2.0mm Small Joint Grasper
242-010-300



2.0mm Small Joint Scissors
242-010-400

Conquest™ Graspers

These instruments are designed to access any area of the shoulder or knee and provide unsurpassed grasping of tissue and loose bodies.



4.0mm Soft Tissue Grasper
300-034-390



3.4mm Claw Rongeur
300-034-380



3.4mm Soft Tissue Grasper
300-034-703



3.4mm Soft Tissue Grasper with Ratchet
300-034-702



3.4mm Grasper
300-034-701


2.7mm Soft Tissue Grasper
300-027-703

2.7mm Soft Tissue Grasper with Ratchet
300-027-702

3.4mm Grasper with Ratchet
300-034-700

Conquest™ Craftsmanship

With increased emphasis on cost containment, Stryker provides a cost-effective instrument with a lifetime warranty on workmanship and materials.



Sterilization Tray (12 and 24 Instruments)




Soft Tissue Grasper with Ratchet

Conquest™ Tip Variations



Straight Tip



15° Up Tip



30° Left Tip



45° Right and Left Tip



Upshaft Tip



30° Right Tip

Conquest™ Features

Low Profile tips allow for precise positioning in extremely tight areas.



All Punches

Stryker's patented Meniscal Crescent design allows the surgeon to take contoured bites of the meniscus.



3.4mm Meniscal Crescent Punches

Stryker's 90° Rotary Punches enable the surgeon to access the anterior to medial third of the meniscus.



3.4mm 90° Rotary Punches

Ordering Information

Punches

2.7mm Straight Big Bite Punch
300-027-100

2.7mm 30° Right Big Bite Punch
300-027-101

2.7mm 30° Left Big Bite Punch
300-027-102

2.7mm 15° Up Big Bite Punch
300-027-103

2.7mm Upshaft Big Bite Punch
300-027-104

2.7mm Narrow Punch
300-027-200

3.4mm Straight Big Bite Punch
300-034-100

3.4mm 30° Right Big Bite Punch
300-034-101

3.4mm 30° Left Big Bite Punch
300-034-102

3.4mm 15° Up Big Bite Punch
300-034-103

3.4mm Upshaft Big Bite Punch
300-034-104

3.4mm 45° Left Tip Biter
300-034-110

3.4 mm 45° Right Tip Biter
300-034-111

3.4mm Narrow Punch
300-034-200

3.4mm 15° Up Narrow Punch
300-034-203

3.4mm Straight Scoop Jaw Punch
300-034-300

3.4mm 30° Right Scoop Jaw Punch
300-034-301

3.4mm 30° Left Scoop Jaw Punch
300-034-302

3.4mm 90° Right Rotary Punch
300-034-400

3.4mm 90° Left Rotary Punch
300-034-401

3.4mm Right Meniscal Crescent Punch
300-034-500

3.4mm Left Meniscal Crescent Punch
300-034-501

3.4mm 30° Left Angled Meniscal Crescent Punch
300-034-502

3.4mm 30° Right Angled Meniscal Crescent Punch
300-034-503

3.4mm Backbiter
300-034-800

Scissors

2.7mm Micro Scissors Punch
300-027-603

3.4mm Hook Scissors
300-034-600

3.4mm 30° Right Hook Scissors
300-034-601

3.4mm 30° Left Hook Scissors
300-034-602

3.4mm Micro Scissors Punch
300-034-603

Graspers

2.7mm Soft Tissue Grasper with Ratchet
300-027-702

2.7mm Soft Tissue Grasper
300-027-703

3.4mm Grasper with Ratchet
300-034-700

3.4mm Grasper
300-034-701

3.4mm Soft Tissue Grasper with Ratchet
300-034-702

3.4mm Soft Tissue Grasper
300-027-703

3.4mm Claw Rongeur
300-034-380

4.0mm Soft Tissue Grasper
300-034-390

Small Joint Instruments

2.0mm Small Joint Punch
242-010-200

2.0mm Small Joint Grasper
242-010-300

2.0mm Small Joint Scissors
242-010-400

Instrument Trays

12 Instrument Sterilization Tray
242-000-012

24 Instrument Sterilization Tray
242-000-024

Joint Replacements

Trauma

Spine

Micro Implants

Orthobiologics

Instruments

Interventional Pain

Navigation

Endoscopy

Communications

Patient Handling Equipment

EMS Equipment

5900 Optical Court
San Jose, CA 95138
t: 800.624.4422 f: 800.729.2917

www.stryker.com

The information presented in this brochure is intended to demonstrate a Stryker product. Always refer to the package insert, product label, and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with [™] designation are trademarks of Stryker.
Products referenced with * designation are registered trademarks of Stryker.

Literature number: 1000-900-260D
JJ/BL 5k 06/04

Copyright © 2004 Stryker
Printed in USA