

**Endoscopy** 

## On The

# **Inside Tract**

FlexVision™ U-500 Flexible Ureteroscope






## **FlexVision**

On the Inside Tract

## **Sophistication**

Superior resolution, deflection, and durability unite in the Stryker Flexible Ureteroscope to provide outstanding visualization and lasting workability. Applying the precision of increased fiber optic bundles, the FlexVision™ U-500 is the market frontrunner in resolution. The latest in flexible optics technology offers surgeons a powerful method to access and visualize the entire kidney, especially the lower pole.

Merging the matchless clarity of the FlexVision™ U-500 with the standardsetting Stryker High Definition camera will provide an unrivaled visual image in the operating room. Urologists will find the Stryker Flexible Ureteroscope an exceptional tool by which to raise the bar in their procedures.



**Matchless Clarity** 

Detailed brightness in the smallest spaces enhances visibility and increases performance during procedures. By joining the forces of the FlexVision™ U-500 with the Stryker Camera System, a surgeon is equipped with the finest tools available to complete complex operations.

## **Intuitive Design**

The lightweight and ergonomic form of the FlexVision™ U-500 makes ureteroscopies simpler and more comfortable. Combining smooth deflection buttons with a sleek handle offers a firm, dependable grip during procedures. During secondary deflection, a locking mechanism allows the hand to rest while accessing the lower pole.

FlexVision™ U-500 Features


- 6.9 Fr distal tip provides non-invasive access.
- 3.6 Fr working channel for optimal treatment and fluid management.
- 275° tight radius deflection.
- Increased fiber optics for high resolution.
- Reinforced proximal end for axial stability and greater surgeon control.

### **Superior Deflection**

Optimal accessibility is provided in both kidneys through primary and secondary deflection to both sides. Kidney treatment is made possible by way of dual deflection with a tight radius, allowing an excellent view of the upper and lower poles.



Dual Deflection Up



Dual Deflection Down

## **Ordering Information**

FlexVision™ U-500 503-888-000

The following products are included:

Cleaning Brush 503-888-003

Lens Cleaner 503-888-004

**Leak Tester** 503-888-005

Port Seal Valve (1 pkg/5) 503-888-001

## **Optional Accessories**

Sterilization/Storage Tray 233-032-800

### Baskets

Stryker also provides a complete line of Disposable Baskets. Please contact your local Stryker Representative for selection and pricing.



Joint Replacements	
Trauma	
Spine	
Micro Implants	
Orthobiologics	
Instruments	
Interventional Pain	
Navigation	
Endoscopy	
Communications	
Patient Handling Equipment	
EMS Equipment	

5900 Optical Court San Jose, CA 95138 t: 800.624.4422 f: 800.729.2917

www.stryker.com

The information presented in this brochure is intended to demonstrate a Stryker product. Always refer to the package insert, product label, and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with  $^\circ$  designation are trademarks of Stryker. Products referenced with  $^\circ$  designation are registered trademarks of Stryker.

Literature number: **1000-900-493C** JJ/BL 5k 01/05