


stryker®  
ENDOSCOPY

# 12K SHAYER



The Power You Want  
The Reliability You Demand






# ULTIMATE BURS


**BARREL BURS-  
Hollow (6 Flute)**  
For rapid, aggressive subacromial decompression, notchplasty and distal clavicle resection without clogging.




**BARREL BURS-  
Hollow (12 Flute)**  
For smoother bone removal



**ROUND BURS-  
Hollow (6 Flute)**  
For rapid, aggressive general abrasion orthoplasty, notchplasty and decompression. New design allows for aggressive bone removal without clogging.



**ROUND BURS-  
Hollow (12 Flute)**  
For smoother bone removal




**EGG BUR-  
Hollow (6 Flute)**  
For rapid, aggressive subacromial decompression, notchplasty and distal clavicle resection without clogging.



**SLAP™ BUR-  
Hollow (6 Flute)**  
The protected bur design allows safe access to the glenoid without compromising the integrity of the labrum.



**UNHOODED  
ROUND BUR**  
For subchondral drilling.



# SMALL JOINT

For use with the Small Joint Handpieces (275-601-000), (275-601-000).

**AGGRESSIVE CUTTER**



**FULL RADIUS**



**HOODED ABRASION**



■ 4.0mm 275-941-000  
■ 5.5mm 275-951-000

■ 3.5mm 275-930-000  
■ 4.0mm 275-940-000  
■ 5.5mm 275-950-000

■ 5.5mm 275-950-012

■ 5.5mm 275-952-000

■ 4.0mm 275-941-500

Developed in conjunction with  
Craig Morgan, M.D.

■ 3.0mm 275-831-000

■ 2.5mm 275-628-000  
■ 3.5mm 275-638-000

■ 2.5mm 275-627-000  
■ 3.5mm 275-637-000

■ 2.0mm 275-641-000  
■ 3.0mm 275-647-000

**Stryker Endoscopy**

2590 Walsh Avenue • Santa Clara, California 95051 • 1.800.624.4422 • Fax 1.800.729.2917

905.332.3235 Canada • 616.385.2600 Middle East • 305.821.1888 Latin America 1000-900-267 Rev. A Made in the U.S.A.



# STANDARD CUTTERS

Hubs color-coded for size only


<p>■ 5.0mm 275-552-000</p>	<p>■ 3.5mm 275-531-000                  ■ 4.0mm 275-541-000                  ■ 5.0mm 275-551-000</p>	<p>■ 3.5mm 275-734-000                  ■ 4.0mm 275-744-000                  ■ 5.0mm 275-754-000                  ■ 5.6mm 275-764-000</p>	<p>■ 4.0mm 280-744-200</p>	<p>■ 4.0mm 275-746-000                  ■ 5.0mm 275-756-000                  ■ 5.6mm 275-766-000</p>	<p>■ 3.5mm 275-732-000                  ■ 4.0mm 275-742-000                  ■ 5.0mm 275-752-000                  ■ 5.6mm 275-762-000</p>	<p>■ 4.0mm 275-748-000                  ■ 4.0mm 275-884-000</p>	<p><b>RESECTOR</b>                  For general debritement. Particularly effective on articular cartilage. Ideal for removal of pterosternum from the acromion. Run in forward for bone removal in acromioplasty and notchplasty.</p>	<p><b>COUGAR</b>                  Open ended design provides both side and end cutting of meniscus, ACL/PCL stump, and plica. Ideal for posterior tears.</p>	<p><b>AGGRESSIVE MENISCUS</b>                  For meniscectomy and articular cartilage debritement.</p>	<p><b>CONVEX ANGLED AGGRESSIVE MENISCUS</b>                  For meniscectomy and articular cartilage debritement. Cutting window on outside.</p>	<p><b>SCALLOPED</b>                  For tough cartilage and ACL stump debritement.</p>	<p><b>FULL RADIUS</b>                  For rapid synovial debritement and meniscus cutting. Excellent for smoothing of articular cartilage.</p>	<p><b>DOUBLE BITE</b>                  Shortens time for extensive synovectomy.</p>	<p><b>END CUTTER</b>                  Open ended for meniscal tear debritement.</p>	<p><b>TRUE END</b>                  For direct end cutting in the posterior compartments.</p>	<p><b>SLOTTED WHISKER</b>                  For safe, selective debritement of chondromalacia and general selective tissue debritement.</p>	<p><b>WHISKER CUTTER</b>                  Smooths tissue edges.</p>	
----------------------------	--------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------	----------------------------	--------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------	-----------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------	-------------------------------------------------------------------------	--



**Stryker Endoscopy**

# Disposable Shaver Blades

## ULTIMATE CUTTERS



### AGGRESSIVE PLUS

Razor sharp toothed inner blade provides aggressive resection of meniscus, cartilage, and synovium. No teeth on outer window to protect articular cartilage upon entry.

- 4.0mm 275-544-000
- 5.5mm 275-564-000



### ANGLED AGGRESSIVE PLUS

New design provides rapid tissue removal without clogging. Unique curved design allows increased access to patella, posterior meniscus, and other hard-to-reach areas. No teeth on outer window to protect articular cartilage upon entry.

- 4.0mm 280-544-100



### TOMCAT CUTTER

The narrow gap between inner and outer cutting surfaces provides aggressive resection of soft tissue in knee & shoulder. Extremely effective on meniscus, & plicae.

- 4.0mm 275-545-000
- 5.5mm 275-565-000



### ANGLED TOMCAT

New design provides rapid tissue removal without clogging. Unique curved design allows increased access to patella, posterior meniscus, and other hard-to-reach areas.

- 4.0mm 280-545-100



### RESECTOR

For general debridement. Particularly effective on articular cartilage. Ideal for removal of periosteum from the acromion. Run in forward for bone removal in acromioplasty and notchplasty.

- 4.0mm 275-542-000
- 5.5mm 275-562-000