


On The **Inside Tract**



Flexibility, Dilation, and Control

Combining a lightweight, ergonomic handle with one-handed thumb control, Stryker Stone Baskets bring stone control and capture to a higher level. Stryker Stone Baskets enhance control and maneuverability with handle-to-basket sheath bonding and strain relief. Combining proven basket configurations with optimum dilation, Stryker is unquestionably the urology solution for maximum control in the operating room.




Nitinol Basket

The Flexibility of Nitinol

Stryker's Nitinol 4 Wire Basket impedes kinking and provides increased distention force. The symmetry, pronounced radiuses, and parachute shape of the basket reduce distortion and trauma to upper tract tissue, while the 1mm round ball Minimum Projecting Tip provides a smooth interface with the working channels of the Stryker flexible scopes. Stryker offers greater strength and tip deflection with the Nitinol Basket's Variable Flexibility Wire Braided Sheath Reinforcement, allowing for consistent and firm wire placement.

The Flat Wire Solution

Stryker's Flat Wire Basket is made of stainless steel flat wire for increased distention force. The proprietary process of shaping and forming baskets facilitates capturing and extracting urinary stones from the ureter and UPJ. Exceptional firmness, effective expansion force, and consistent wire positioning are all within the surgeon's grasp with Stryker's symmetrical 4 Flat Wire basket design.



Flat Wire Basket

Ordering Information

Flat Wire 4 Wire Baskets

Order Number	Size	Length	Basket Size	Tip
504-824-090	2.4Fr	90cm	16mm	Round
504-824-115	2.4Fr	115cm	16mm	Round
504-803-070	3.0Fr	70cm	10mm	Round
504-803-090	3.0Fr	90cm	20mm	Round
504-803-115	3.0Fr	115cm	20mm	Round

3 Wire Grasping Forceps Baskets

Order Number	Size	Length
504-833-070	3.0Fr	70cm
504-833-090	3.0Fr	90cm
504-833-115	3.0Fr	115cm

Nitinol 4 Wire Baskets

Order Number	Size	Length	Basket Size	Tip
504-881-070	1.9Fr	70cm	10mm	MPT
504-881-090	1.9Fr	90cm	10mm	MPT
504-881-115	1.9Fr	115cm	10mm	MPT
504-882-090	2.4Fr	90cm	12mm	MPT
504-882-115	2.4Fr	115cm	12mm	MPT
504-883-070	3.0Fr	70cm	10mm	MPT
504-883-090	3.0Fr	90cm	16mm	MPT
504-883-115	3.0Fr	115cm	16mm	MPT

Products referenced with ™ designation are trademarks of Stryker.
Products referenced with ® designation are registered trademarks of Stryker.

Literature number: 1000-900-524A
JJ/BL 5k 7/05

Copyright © 2005 Stryker
Printed in USA

stryker®

5900 Optical Court
San Jose, CA 95138
t: 800.624.4422 f: 800.729.2917

www.stryker.com