

On The Inside Tract

Visual Urethrotome



Visual Urethretome

On the Inside Tract

The Urology Solution

Stryker is a leader in visualization solutions for minimally invasive procedures. EndoUrologists may now combine the leading 3-chip camera and scope technology with durable, easy to use rigid hardware and instrumentation. Efficient procedures and clear medical images are now packaged into one Urology solution, offered by one company: Stryker.

Direct Vision

The Stryker Visual Urethretome is utilized for incision of urethral strictures and bladder neck contractures under direct vision. The Visual Obturator converts the Visual Urethretome to a direct vision urethroscope.

Clear and Crisp

The goal of Stryker autoclavable endoscopes is simple: provide a clear, crisp image for Urologists. Combined with Stryker's new high definition (HD) camera technology, an excellent picture is now accessible through an advanced and comfortable tool. First-rate visualization is ensured through Stryker Rigid Optical instrumentation, a line of hardware created for durability and endurance.

Features

- » 5Fr port accommodates guidewire or small catheter.
- » Visual Obturator for less traumatic sheath insertion.
- » SpeedLock™ locking mechanism assembles quickly and securely.



Visual Obturator



Standard Obturator

Visual Urethretome

21Fr Visual Urethretome Sheath
502-880-521


25Fr Visual Urethretome Outer Sheath
502-880-525

21Fr Visual Urethretome Optical Bridge
502-880-524

21Fr Urethretome Standard Obturator
502-880-522

21Fr Urethretome Visual Obturator
502-880-523

25Fr Balloon Catheter Sheath
502-880-526



Optical Bridge

Accessories


Cold Knife, Half Moon
502-880-501

Cold Knife, Straight
502-880-502

Cold Knife, Serrated
502-880-503

Cold Knife, Hook
502-880-504

Cold Knife, Loop
502-880-505



Cold Knives

Joint Replacements

Trauma

Spine

Micro Implants

Orthobiologics

Instruments

Interventional Pain

Navigation

Endoscopy

Communications

Imaging

Patient Handling Equipment

EMS Equipment

5900 Optical Court
San Jose, CA 95138
t: 800.624.4422 f: 800.729.2917

www.stryker.com

The information presented in this brochure is intended to demonstrate a Stryker product. Always refer to the package insert, product label, and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with ™ designation are trademarks of Stryker.
Products referenced with ® designation are registered trademarks of Stryker.

Literature number: 1000-900-523A
JJ/BL 5k 4/05

Copyright © 2005 Stryker
Printed in USA